



Directions for use

Manufacturer / pharmaceutical company:

Tetra GmbH • Herrenteich 78 • D-49324 Melle • Tel. 05422/105-0 • Made in Germany • www.tetra.net

This veterinary medicine is marketed in accordance with the Small Animal Exemption Scheme.
Manufacturing Authorisation Number SAM 0005.

Tetra Medica FungiStop

Method of administration: Solution for use in aquarium water. 100 ml per original bottle

Amount of active ingredient per 100 ml of solution: colloidal silver 0.240g

Areas of application: Fungal infections caused by Achlya, Saprolegnia, etc., spawn fungus, bacterial mucous membrane infections, mouth decay and mouth rot, infectious wounds, heals wounds safely by providing overall mucous membrane protection. **Disease symptoms:** Cotton wool-like, white patches, mould-like structures on the fish skin, fungal wounds, skin inflammations.

Contraindications and side effects: None known. The owner is requested to inform a veterinary surgeon or pharmacist of any newly-discovered side effect. The specialist retailer or product manufacturer should also be informed of the side effects.

Remedy for freshwater tropical fish

Dosing instructions / method of application: Dispense the product just once on the first day of the treatment cycle, using the scale lines (each line representing 10ml) on the side of the bottle.

For established aquaria:

Normal dosage: 5 ml Tetra Medica FungiStop to 20 l aquarium water
In cases of severe infection: 5 ml to 10 l

For quarantine tanks: 5 ml to 20 l (without using biologically-active filters and equipment)

Shake the product well before use! Do not use any water conditioners that neutralise heavy metals prior to or during treatment. The treatment cycle lasts for 5 days and begins as soon as the one-off dosage is administered on the first day. At the end of each treatment cycle, carry out a partial water change and use activated carbon to filter the water for a period of 24 hours. In the event that no improvement is seen, a second dosage can be added to the tank once the water has been filtered using activated carbon. The product should be added in the normal dosage to prevent spawning of fungus.

Instructions for correct use:

During treatment, please ensure that the aquarium is well aerated and do not use absorbent filter media (such as activated carbon, synthetic resin, etc.) to filter the water. Switch off the UV steriliser. Tetra Medica FungiStop will not damage plant life. Do not use Tetra Medica FungiStop in conjunction with any other remedies.

Special storage precautions:

Keep the product away from sunlight! Do not store at temperatures above 25°C! The shelf life after the original packaging has been opened is 6 months.

Special precautions:

For Animal Treatment only. Keep this remedy out of the reach and sight of children! Do not use the remedy once the expiry date has passed! Do not use the product in tanks containing invertebrates (e.g. snails, shrimps and other crustaceans)!

Dispose of this product and its container at a hazardous or special waste collection point.

Information produced: August 2011